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
**SHUTTING VALVES
VANNES D'ENTRANGLEMENT
DROSSELVENTILE
VALVOLAS PARCIALIZADORAS**

- MANUALE ISTRUZIONI
- OPERATOR'S HANDBOOK
- MANUEL D'UTILISATION
- GEBRAUCHSANWEISUNG
- MANUAL DE INSTRUCCIONES

CATTANI S.p.A.

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INDICE

Pagina

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INDEX

Page

— FUNCTIONING GENERAL DATA	8
— INTRODUCTION	8
— SIGNS AND WARNINGS	8
— ELECTROPNEUMATIC VALVES	9
— FUNCTIONING	9
— INSTALLATION	9
— GENERAL RECOMMENDATIONS AND BIOLOGICAL HAZARD	9
— MAINTENANCE	10
— IMPORTANT NOTES	10
— TRANSPORT AND STORAGE	10
— TRANSPORT OF SECOND-HAND APPLIANCES	11

Functioning general data

Electropneumatic valve AC/DC 24 V Mignon electropneumatic valve 04

	Electropneumatic Valve AC/DC 24 V	Electropneumatic Valve Mignon 04
Rated voltage:	24 V $\cong$ AC/DC	24 V $\cong$ AC/DC
Rated frequency:	50/60 Hz	50/60 Hz
Rated current:	0,33 A $\cong$ AC/DC	0,12 A $\cong$ AC/DC
Used power:	8 W	3 W
Insulation Class:	I	I
Type of apparatus:	B	B
Use:	continuous working	continuous working
Level of protection against liquids:	common	common
Level of protection against direct or indirect contacts type:	B 	B 
Suction ready signal relay contact features:	Max 5 A 24 V AC Cos. $\varphi = 1$	Max 1 A 24 V AC Cos. $\varphi = 1$

Introduction

The following presentation serves to illustrate the equipment and systems to professionals and repair experts. **We think it will also be useful to explain their operation and maintenance and inform both experts and users about the hazards and precautions that should be taken to prevent accidents.**

Signs and warnings

- **Biological danger, danger of infections from epidemic diseases.**

- **High temperature.**

- **General danger sign.**


- **Compulsory direction of flow and rotation.**


Signs cannot always fully express danger warnings, therefore it is necessary that the user reads the warnings and keeps them in due consideration. Failure to respect a danger sign or warning may harm operator or damage the equipment. Do not remove protections, do not tamper with machines or their operation.

Electropneumatic valves, item 021144 and item 024150 (draw. page 30)

They have been designed to open and close the piping connecting every dental unit serviced by one aspiration system. They work in presence of air, as when they are installed on Mini-Separators (Fig. 1), as well as in presence of aspirated fluid, as with Turbo-Jet modular and Separator Tanks (Fig. 2). A more accurate maintenance is required in this case.

(021144)
**Electropneumatic
Valve**
AC/DC 24 V

(024150)
**Mignon
Electropneumatic
Valve 04**

Functioning (draw. page 31)

With the vacuum function operating, the depression draws the membrane onto the vacuum hose and thereby closes the air passage (part A fig. 3); by raising a terminal from its housing, the solenoid valve is energized and opens a passage that causes a depression on the other side of the membrane as well.

The two depressions are equal and thus the air flow opens the suction (part B fig. 3).

The membrane (Fig. 3 det. C) can work up to a max. head of 3500 mm H₂O.

Installation (draw. page 30-33)

These valves can be installed inside the water unit, on the Mini-Separators (Fig. 1) or, in absence of this one, on the aspiration piping (Fig. 2), at a min. height of 200 mm above the draining pipe running on the floor or grooved under it (Fig. 6); it is recommended to avoid any sharp bends (see dashed line), which make the aspirated fluid rebound.

The valve must be steadily fixed in upright position and fed at 24V.

Aspirated fluid usually passes through the tip support filter; failing this one, it is recommended to use the filter installed on the electrovalve itself (Fig. 2 det. A).

General recommendations and biologic hazard

- **Biological danger, danger of infections from epidemic diseases.**

- **General danger sign.**

Before servicing any used equipment, carry out some washing operations with Puli-Jet (Fig. page 10), or other disinfectant appropriate for this purpose. **Disconnect electric supply** and padlock mains switch if provided.

Before servicing the equipment wear disposable gloves, goggles, mask.

Aspirated debris is infected and contaminated, moreover the appliances to be serviced could be under pressure and the risk of contaminating splashes is clearly higher.

Maintenance (draw. page 31-32)

- **Biological danger, danger of infections from epidemic diseases.**

- **General danger sign.**

After every working day clean the system by aspirating Puli-Jet, diluted in hot water (50 °C) according to instructions. Aspirate the solution by means of Pulse Cleaner, which will add the necessary turbulence to the detergent effect of Puli-Jet, to optimise the cleaning (see side illustration).


Use of products different from those recommended by the manufacturer can jeopardise the functioning of the membrane (Fig. 3, det. C) and cause damages.

Clean filters on tip support and/or on valve (Fig. 2, det. A).

Technical inspections should be made quarterly by an expert authorized by Cattani S.p.A. and in possession of original replacement parts.

Before operating on the valve check the exploded drawing (fig. 5). Check the inside of the valve for cleaning (inspect channel, part D, fig. 4) and if necessary remind the clinic personnel to use Puli-Jet and Pulse Cleaner.

By energizing the solenoid valve it is possible to control the movement of the magnet, the membrane and the relays involved. Replace any misshapen membrane.

Puli-Jet and Pulse Cleaner

Important notes

- The manufacturer is at your disposal to supply spare parts, literature, instructions and anything that might be useful.
- Our Concessionaires, Agents, Retailers, and authorised Service Engineers are provided with split up drawings, electrical diagrams, instructions and updating about maintenance and service.
- The guarantee and manufacturer's liability cease in case the appliances and/or the plants have been tampered with operations of any kind done by persons unfit and so unauthorised by the manufacturer.
- For any usage not contemplated or specified in this manual, please contact the manufacturer.

Transport and storage

- Packed equipment can be transported and stored at a temperature range of - 10 °C + 60 °C.
- Packages must be kept away from water and splashing and cannot tolerate humidity >70%.

- *Packages with the same weight can be stored in piles of three only.*

Transport of second-hand appliances

- *Prior to packing, clean and sanitize with Puli-Jet (see “General recommendations and biologic hazard” and “Maintenance”).*
- *Place the appliance into a polyethylene bag, seal it and pack in a 3-layer corrugated board.*

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ESAM S.p.A.